


# S-2765™

## CHROMOGENIX

Instrumentation Laboratory Company - Bedford,  
MA 01730-2443 (USA)

Instrumentation Laboratory SpA-V.le Monza 338-  
20128 Milano (Italy)



ENGLISH

# S-2765™

For Laboratory Use Only  
For General Laboratory Use

S-2765 is a chromogenic substrate for determination of Factor Xa. It is also very sensitive to trypsin.

#### COMPOSITION

Each vial contains the chromogenic substrate S-2765, 25 mg and mannitol 60 mg as a bulking agent.

#### CHEMISTRY

**Chemical name:** N- $\alpha$ -Benzoyloxycarbonyl-D-arginyl-L-glycyl-L-arginine-p-nitroaniline-dihydrochloride

**Formula:** N- $\alpha$ -Z-D-Arg-Gly-Arg-pNA · 2HCl

**Mol. wt.:** 714.6

**$\epsilon_{316\text{ nm}^2}$ :**  $1.27 \cdot 10^4 \text{ mol}^{-1} \cdot \text{L} \cdot \text{cm}^{-1}$

**Solubility:** > 40 mmol/L in H<sub>2</sub>O  
> 10 mmol/L in Tris buffer (pH 8.3, I 0.25)

**Stability:** Lyophilized substance: stable at 2-8°C until expiry date printed on the label. The substance is hygroscopic and should be stored in a dry place.  
Solution: 2 mmol/L in H<sub>2</sub>O is stable for six months at 2 to 8°C  
Contamination by microorganisms may cause hydrolysis.

**Suitable stock solution:** 1-2 mmol/L in H<sub>2</sub>O.

#### PRINCIPLE

$\text{N-}\alpha\text{-Z-D-Arg-Gly-Arg-pNA} \xrightarrow{\text{Enzyme}} \text{N-}\alpha\text{-Z-D-Arg-Gly-Arg-OH} + \text{pNA}$

The method for the determination of activity is based on the difference in absorbance between the pNA formed and the original substrate. The rate of pNA formation, i.e. the increase in absorbance per second at 405 nm, is proportional to the enzymatic activity and is conveniently determined with a photometer.

#### KINETIC DATA

**Factor Xa (bovine):**  $k_m = 1 \cdot 10^{-4} \text{ mol/L}$ ,  $k_{cat} = 290 \text{ sec}^{-1}$  in Tris buffer pH 8.3, I 0.25 at 37°C.

**Factor Xa:** (human plasma activated with Russel's Viper Venom):  
 $k_m = 3 \cdot 10^{-4} \text{ mol/L}$  in Tris buffer pH 7.8, I 0.4 at 37°C.

#### STANDARDIZATION

An activity of  $\Delta A/\text{min} = 0.05$  (37°C) is obtained by using a substrate concentration of  $2 \cdot k_m$  and a concentration of 0.04 nkat/mL of FXa (Chromogenix AB).



## CHROMOGENIX