


Tissue Plasminogen Activator (t-PA)

Method Sheet S-2288

Determination of Tissue Plasminogen Activator (t-PA) in Purified Preparations

Principle:



The tissue plasminogen activator (t-PA) is a family of serine proteases which activate plasminogen by splitting a single Arg-Val bond of the plasminogen molecule. In purified systems these enzymes have been shown to hydrolyse tripeptide chromogenic substrates (1, 2, 3).

The t-PA activity is thus determined by the rate at which p-nitroaniline (pNA) is released. The formation of pNA can be followed spectrophotometrically at 405 nm by using a recorder (initial rate method).

The correlation between the change in absorbance per minute ($\Delta A/\text{min}$) and the t-PA activity is linear in the 0.05 - 0.5 $\mu\text{kat/l}$ or 3 - 30 U/l range.

The amidolytic activity does not necessarily parallel the fibrinolytic activity for different t-PA preparations.

Reagents:

1. Substrate S-2288 (H-D-Ile-Pro-Arg-pNA • 2HCl) mol.wt. 577.6.

The substance (25 mg) is dissolved in 8.65 ml (t-PA one-chain) or 43 ml (t-PA two-chain) of distilled water. If prepared with sterile water and not contaminated by microorganisms, the 5 mmol/l or 1 mmol/l solution is stable for at least two months at 2 to 8°C.

2. Buffer pH 8.4 (25°C)

Tris	12.1 g (100 mmol/l)
NaCl	6.2g (106 mmol/l)
Distilled water	800 ml

Adjust the pH to 8.4 at 25°C by adding an appropriate amount (approximately 44 ml) of 1 mol/l HCl. Fill up to 1000 ml with distilled water. The buffer, if not contaminated, is stable for six months at 2 to 8°C.

3. Acetic acid 20%

Acetic acid is used in the acid-stopped method.

Equipment :

1. Spectro- or filter photometer, 405 nm with cuvette housing, thermostated at 37°C
2. Semi-microcuvettes, 1 cm
3. Thermostat, 37°C
4. Stop watch
5. Disposable plastic tubes

Procedure:

Sample:

Purified tissue plasminogen activator is dissolved in buffer to an enzyme activity of 0.05 - 0.5 $\mu\text{kat/l}$ (3 - 30 U/l). See note*.

It has been advised to use a surfactant to avoid adsorption to surfaces. A final concentration of 0.1 g/l of Triton X-100 (Merck) is recommended.

Method:

Initial rate method		Volume
Buffer	μl	200
Thermostate at 37°C (2-4 min)		X
Sample (20-25°C)	μl	200
Mix and thermostate at 37°C (2-4 min)		X
S-2288 (37°C)	μl	200
Mix		X

Transfer sample immediately to a 1 cm semi-microcuvette (preheated to 37°C) for measurement of the absorbance change in a photometer at 405 nm and at 37°C. Calculate $\Delta A/\text{min}$.

Calculation:

The t-PA activity in the prepared tissue plasminogen activator solution is calculated from the following formulas:

$$\mu\text{kat/l} = \Delta A/\text{min} \times 5.21$$

$$\text{U/l} = \Delta A/\text{min} \times 313$$

***Note:**

In the test (600 μl) 0.25 μg (\cong 100 IU) of the procine heart tissue activator gives:

$$\Delta A/\text{min} \cong 0.012 \text{ (one-chain)}$$

$$\Delta A/\text{min} \cong 0.065 \text{ (two-chain)}$$

Bibliography

Friberger P et al. Activity of plasminogen activators on tripeptide chromogenic substrates. In: Progress in Chemical Fibrinolysis and Thrombolysis. Vol.4. Davidson J.F. et al (eds): Churchill Livingstone, 149-153 (1979).

Friberger P: Chromogenic Peptide Substrates. Scand J Clin Lab Invest 42, suppl. 162 p 57 (1982).

Verheijen JH et al. Quantitative analysis of the composition of mixtures of one-chain and two-chain tissue-type plasminogen activator with a spectrophotometric method. Thromb Res 39, 281-288 (1985).

Dodd I et al. Isolation and preliminary characterisation of active B-chain of recombinant tissue-type plasminogen activator. Thromb Haemost 55, 94-97 (1986).

Rijken DC and Groenveld E. Isolation and functional characterisation of the heavy and light chains of human tissue-type plasminogen activator. J Biol Chem 261, 3098-3102 (1986).

Mitsubayashy S et al. Plasminogen activator in bladder tumors. Urol Res 15, 335-339 (1987).

Sugawara Y et al. Kinetic analysis of the enhancement of the activities of t-PA induced by the presence of monoclonal antibody (C9-5). Thromb Res 50, 637-646 (1988).

Urano T et al. Stimulation of the amidolytic activity of single chain tissue-type plasminogen activator by fibrinogen degradation products: possible fibrin binding sites on single chain tissue-type plasminogen activator molecule. Biochim Biophys Acta 1077, 245-252 (1991). Ball EL et al. Selection of monoclonal antibodies that bind and inhibit tissue-type plasminogen activator. Hybridoma 12, 317-326 (1993).

Feng YH and Hart G. In vitro oxidative damage to tissue-type plasminogen activator: a selective modification of the biological functions. Cardiovasc Res 30, 255-261 (1995). Higazi AA et al. Identification of an inhibitor of tissue-type plasminogen activator-mediated fibrinolysis in human neutrophils. A role for defensin. J Biol Chem 270, 9472-9477 (1995).