


Standardisation of Urokinase

Method Sheet

Principle



The urokinase activity is determined by its amidolytic effect on the substrate pyro-Glu-Gly-Arg-pNA. The rate at which p-nitroaniline (pNA) is released is measured photometrically at 405 nm. This can be followed on a recorder (initial rate method) or read after stopping the reaction with acetic acid (acid stopped method).

The correlation between $\Delta A/\text{min}$ (or absorbance) and the urokinase activity is linear in the range 5-40 Ploug or CTA units.

The urokinase concentration should preferably be given in units of substrate splitting activity, but may also be calculated by using standards prepared from a standard urokinase preparation. The amidolytic activity, however, does not necessarily parallel the fibrinolytic activity for different urokinases.

Equipment

- | | | |
|---|-----------------------------|--|
| 1. Spectro- or filterphotometer, 405 nm | 4. Stop-watch | <i>Additional equipment for the initial rate method:</i> |
| 2. Semi-microcuvettes, 1 cm | 5. Disposable plastic tubes | |
| 3. Waterbath or thermostate, 37°C | | 6. Photometer with cuvette housing, thermostated at 37°C |
| | | 7. Recorder |

Reagents

- Substrate S-2444 (pyro-Glu-Gly-Arg-pNA•HCl), mol. wt. 498.9 Art. No. 82 03 57**
Each vial contains 25 mg S-2444. The substance is dissolved in 16.7 ml of distilled water. If prepared with sterile water and not contaminated by micro-organisms, the 3 mmol/l solution is stable for at least six months at 2-8°C.

- Urokinase standard**
The urokinase standard is dissolved in or diluted with Solvent (Reagent 3) to a concentration of 400 units/ml (Ploug or CTA units). The dilution is stable for one day at 2-8°C.

- Solvent**
Distilled water containing 5 g/l of Carbowax 6000 (Union Carbide, N.Y., USA).

- Buffer pH 8.8 (25°C)**

Tris	6.1 g (50 mmol/l)
NaCl	2.2 g (38 mmol/l)
Distilled water	800 ml

Adjust the pH to 8.8 at 25°C by adding an appropriate amount (approximately 12 ml) of 1 mol/l HCl. Fill up to 1000 ml with distilled water. The buffer, if not contaminated, is stable for two months at 2-8°C.

Note:

Although the substrate is quite selective, there may be a risk for influence of other proteases if the preparation is heavily contaminated. The addition of Trasylol® (aprotinin) (Bayer, FRG), 10 KIU/ml, to the buffer may in such cases be favourable.

- Acetic acid 20%**
Acetic acid is used in the acid-stopped method.

Procedure

Sample

The urokinase is dissolved in or diluted with Solvent (Reagent 3) to a concentration of approximately 400 units/ml (Ploug or CTA units).

By using commercially available urokinase (Leo or Abbott) it was found that the dilution was stable for at least one day when kept at 2-8°C.

Note:

If the urokinase preparation is contaminated with proteolytic enzymes, Trasylol® (aprotinin) (Bayer, FRG) may be added to a concentration of 10 KIU/ml in order to increase the stability.

Standardisation

40 units: Use the urokinase standard 400 units/ml (Reagent 2). (See Methods)

5 units: Use the urokinase standard 400 units/ml (Reagent 2) diluted 1:8 with buffer (Reagent 4). (See Methods)

Standard curve: See appendix.

Methods

Initial rate method		Sample
Buffer	μl	800
Thermostate at 37°C (5-10 min)		X
Urokinase sample or urokinase standards	μl	100
Mix and thermostate at 37°C (1-2 min)		X
S-2444 (37°C)	μl	100
Mix		X

Transfer the sample immediately to a 1 cm semi-microcuvette (preheated to 37°C) for measurement of the absorbance change in a photometer at 405 nm and at 37°C. Calculate $\Delta A/\text{min}$.

Acid stopped method		Sample	Sample blank
Buffer	μl	800	800
Thermostate at 37°C (5-10 min)		X	-
Urokinase sample or urokinase standards	μl	100	100
Mix and thermostate at 37°C (1-2 min)		X	-
S-2444 (37°C)	μl	100	-
Mix and incubate		300	-
Acetic acid 20%	μl	100	100
Mix immediately		X	X
S-2444 (37°C)	μl	-	100
Mix		-	X

Read the absorbance (A) of the sample against a water or sample blank in a photometer at 405 nm. The colour is stable for at least 4 hours.

Calculation

Plot $\Delta A/\text{min}$ or A for the standards against their known urokinase activity.

Calculate the urokinase activity of the sample in Ploug or CTA units. By multiplying the results with 10 the concentration in units/ml is obtained.

The urokinase activity can also be calculated from the following formulas:

$$\begin{aligned}\text{Initial rate method: } \mu\text{kat/l} &= \Delta A/\text{min} \times 17.4 \\ \text{U/l} &= \Delta A/\text{min} \times 1042\end{aligned}$$

$$\begin{aligned}\text{Acid stopped method: } \mu\text{kat/l} &= A \times 3.8 \\ \text{U/l} &= A \times 229\end{aligned}$$

References

1. Claeson G et al. Methods for determination of prekallikrein in plasma, glandular kallikrein and urokinase. *Haemostasis* 7, 76-78 (1978).
2. Paar D and Marhuln D. Spectrophotometric determination of urokinase in urine after gel filtration, using the chromogenic substrate S-2444. *J Clin Chem Clin Biochem* 18, 557-562 (1980).
3. Friberger P: Chromogenic Peptide Substrates. *Scand J Clin Lab Invest* 42, suppl. 162, 55 (1982).
4. Philo R D and Gaffney P J: Relative potencies of different molecular weight forms of urokinase. In: *Progress in Chemical Fibrinolysis and Thrombolysis*. Vol. V Ed Davidson J F, Nilsson I M and Astedt B, 220-222 (1980).
5. Millar WT and Smith JF. Comparative study of different urokinase preparations. *Thromb Res* 30, 425-429 (1983).
6. Higazi AA and Mayer M. In vitro inhibition of urokinase by penicillins. *Thromb Haemost* 60, 305-307 (1988).
7. Keabian PR and Henkin J. A chromogenic enzymatic assay capable of detecting prourokinase-like activity material in plasma. *Thromb Res* 65, 401-407 (1992).
8. Bhargava K and Ando HY. Immobilization of active urokinase on albumin microspheres: use of a chemical dehydrant and process monitoring. *Pharm Res* 9, 776-781 (1992).

Appendix

Standard curve

The urokinase standard 400 units/ml (Reagent 2) is further diluted according to the table below.

Ploug or CTA units (see Methods)	Urokinase standard 400 units/ml μ l	Solvent μ l
5	100	700
10	100	300
20	200	200
30	300	100
40	400	-