BIOPHEN™ CS-11(22)

Factor Xa Chromogenic Substrate

Vial of 25 mg/flacon de 25 mg

Ref. 229015

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Synthetic generic Chromogenic substrate, highly purified and stabilized.

Specificity:

Recommended substrate for Factor Xa.

Presentation:

Vial containing 25 mg of BIOPHEN™ CS-11(22), lyophilised in presence of Mannitol as a bulking agent.

Reconstitution:

According to the research protocol used, the BIOPHEN™ CS-11(22) chromogenic substrate can be restored with variable volumes of distilled water ; for example 10 mL can be used for a substrate concentration of 2.5 mg/mL, or 20 mL for a substrate concentration of 1.25 mg/mL. Shake thoroughly until complete dissolution (vortex). Let to stabilize for 30 min. at room temperature.

Peptide sequence:Mixture (50%-50%) of Bz-Ile-Glu (γ OCH₃)-Gly-Arg-pNa (form 1) and Bz-Ile-Glu (γ OH)Gly-Arg-pNa (form 2)**Purity grade:**

> 95%.

Molecular weight:

711.8 (form 1) and 697.7 (form 2).

pNA content (measured on RP-HPLC):

RP-HPLC < 0.1%

Use:All research studies and protocols where a source of chromogenic substrate for Factor Xa is required. **This kit is for research use only and should not be used for patient diagnosis or treatment.****Respective reactivities:**

FXa	Thrombin	Plasmin	Kallicrein	aPC
100	1	2	3	0

Assay conditions must be duly established for rendering the assay conditions totally specific for Factor Xa, when this substrate is used.

Stability:

- Restored :**
- 7 days at room temperature
 - 3 months at 2-8 °C
 - **Do not freeze.**

Lyophilized : - Until the expiration date printed on the packaging.

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