# Anti-Thrombin (AT)

REF PP004B, Pn 1 x 1 mL

REF PP004L, Pn 10 x 1 mL

Anti-Thrombin (AT)  
**FOR RESEARCH USE ONLY.**  
**NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

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English, last revision 07-2020

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**INTENDED USE:**

All research studies and protocols where a source of highly purified human AT is necessary.  
**This kit is for research use only and should not be used for patient diagnosis or treatment.**

**REAGENTS:**

**[Pn]** Highly purified from human plasma using a combination of affinity and ion-exchange chromatographies.

Vial containing at least 9 IU of AT (when determined and standardized against the WHO/NIBSC International standard for Antithrombin, Concentrate, human).

**1 vial of 1 mL (PP004B).**  
**10 vials of 1 mL (PP004L).**

- Excipients:** Content per vial: Glycine (20 mg), Hepes (6 mg), Sodium Chloride (9 mg).
- Purity:** One major band at about 58,000 daltons on SDS-PAGE.
- Activity:** In a clotting assay for the thrombin inhibition capacity in presence of an excess of unfractionated heparin : 9 IU AT neutralize at least 2,000 NIH units of thrombin. Tested for anti-Xa activity in presence of heparin, with a chromogenic assay (> 6 IU/mg). Factor Xa neutralisation capacity is > 500 µg of Xa/mg AT. The absence of heparin is verified after purification, using a toluidine blue assay (refer to: F.C MacIntosh, "A colorimetric method for the standardization of heparin preparations", Biochem, 1941, 35:776-782), and an anti-Xa specific chromogenic assay.

**WARNING AND PRECAUTIONS:**

- Some reagents provided in these kits contain materials of human origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- Waste should be disposed of in accordance with applicable local regulations.
- Aging studies show that the reagents can be shipped at room temperature without degradation.
- This device of *in vitro* use is intended for professional use in the laboratory.

**REAGENT PREPARATION:**

Gently remove the freeze-drying stopper, to avoid any product loss when opening the vial.

**[Pn]** Reconstitute the contents of each vial with exactly **1 mL of distilled water** in order to get a concentration of about 9 IU/mL. Dilute in the adequate buffer, if necessary. Shake vigorously until fully dissolved. Homogenize and allow to stabilize for 15 minutes at room temperature (18-25°C), homogenize before use.

**STORAGE AND STABILITY:**

Unopened reagents should be stored at 2-8°C in their original packaging. Under these conditions, they can be used until the expiry date printed on the kit.

**[Pn]** Reagent stability after reconstitution, free from any contamination or evaporation, and stored closed, is of:

- 1 month** at 2-8°C.
- 72 hours** at room temperature (18-25°C).
- 6 months** frozen at -20°C or less\*

\*Thaw only once, as rapidly as possible at 37°C, adapting the incubation period to the volume of reagent. The stability of the thawed reagent should be checked under laboratory work conditions. (tested for anti-Xa chromogenic activity in presence of heparin: > 6 IU/mg).

**LIMITATIONS:**

- Any reagent presenting an unusual appearance or showing signs of contamination must be rejected.

**The results obtained should be for research purposes only and not used for patient diagnosis or treatment.**

**SYMBOLS:**

Symbols used and signs listed in the ISO 15223-1 standard, see Symbol definitions document.

*Changes compared to the previous version.*

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