Bovine Thrombin

Thrombine bovine

21 NIH/ Vial / 21 NIH/flacon (BE102A)
Set of 6 vials / Coffret de 6 flacons (BE102K)

Ref. BE102A/BE102K

STORE AT
2-8°C

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC
PROCEDURES.

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English, last revision 03-2025

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ORIGIN:

Prepared from a prothrombin complex concentrate extracted from bovine citrated plasma, and fully activated. Thrombin is then purified using an ion exchange chromatography and stabilised. It is a highly purified preparation, mainly in the α form, with a specific activity ≥ 900 NIH/mg.

PRESENTATION:

Thrombin activity is reported in harmonized NIH or IU or USP units by reference to the WHO/NIBSC International Standard for Human Thrombin (19/188) (refer to "A reunification of the US ("NIH") and International Unit into a single standard for Thrombin", Longstaff et al, Thromb Haemost 2005, 93:261-6)). NIH is a clotting unit for thrombin activity.
Stabilised and lyophilised bovine Thrombin.
Vial containing approximately 21 NIH of thrombin.
BE102A: 1x21 NIH; BE102K: 6x21 NIH.
This kit is for research use only and should not be used for patient diagnosis or treatment.

RECONSTITUTION:

Each vial must be restored with 1 mL distilled water, in order to get a concentration of about 21 NIH/mL. Dilute in the adequate buffer, if necessary.

EXCIPIENTS:

Content per vial: Poly-Ethylene-Glycol 6000 (PEG 6000) (20mg), Sodium Chloride (9mg), Glycine (20mg), BSA (10mg), stabilizers.

PURITY:

Bovine Thrombin has one major band of 35,000 daltons on SDS-PAGE.

ACTIVITY:

Clotting time after reconstitution with 1mL and dilution at 1:3:
on purified bovine fibrinogen (5 mg/mL) : ≤ 12 sec"
on purified human fibrinogen (5 mg/mL) : ≤ 12 sec"
Chromogenic activity on CS-01(38) thrombin substrate: A405 : > 1.20 /min.NIH.mL
Clotting activity determined in NIH (or IU).
Chromogenic activity on specific substrates (expressed in nkats): tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the specific activities are defined on the various IIa substrates (expressed in nkats/ μ g):

CS-01(38)	CS-01(81)
-----------	-----------

The exact activity in nkats/ μ g with CS-01(38) is reported on the analysis certificate for each lot.

VIRAL SAFETY:

The bovine plasma used for Bovine Thrombin preparation was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

USE:

For *in vitro* use only. All research studies and protocols where a source of highly purified bovine Thrombin is necessary.

The results obtained should be for research purposes only and not used for patient diagnosis or treatment.

STABILITY:

Restored : - 2 days at room temperature (18-25°C).
- 7 days at 2-8°C.
Lyophilized : Until the expiration date printed on the kit.

The stability of the reconstituted reagent should be checked under laboratory work conditions.

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