# tPA Stimulator

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Vial of 1 mg/Flacon de 1 mg (Ref PP021B)  
Set of 6 vials of 1 mg/ Coffret de 6 flacons de 1 mg (Ref PP021K)

Ref. PP021B / PP021K

STORE AT  
2-8 °C

For research use only

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CONSERVER A  
2-8 °C

### Origin:

Tissue Plasminogen Activator (tPA) stimulator, lyophilised reagent obtained with oxidized fibrin(ogen) fragments according to the method of Stief et al (1).

### Presentation:

Vial containing approximately 1 mg of tPA Stimulator.

### Reconstitution:

Each vial must be restored with 2 ml distilled water. It can be diluted to the desired concentration when required. Use a buffer containing a carrier protein (such as BSA) for dilution, in order to avoid material loss.

### Excipients:

- Phosphate, Mannitol.
- About 17.5 mg sodium chloride/vial

### Viral safety:

Fibrinogen used for preparing the tPA stimulator was extracted from human plasma, tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

### Use:

For *in vitro* use only.  
All research studies and protocols where a source of tPA Stimulator is necessary.

### Stability:

- Restored :
- 72 hours at room temperature (18-25 °C)
  - 7 days at 2-8 °C
  - 6 months frozen at -20 °C or below

Lyophilized : Until the expiration date printed on the vial.

Reference (1): Stief et al, "trombin converts singlet oxygen (1O<sub>2</sub>) -oxidized Fibrinogen in a soluble tPA cofactor. A new method for preparing a stimulator for functional tPA assays"; Ann Hematol, 2001, 80:189-194.

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Reference (1): Stief et al, "trombin converts singlet oxygen (1O<sub>2</sub>) -oxidized Fibrinogen in a soluble tPA cofactor. A new method for preparing a stimulator for functional tPA assays"; Ann Hematol, 2001, 80:189-194.

# tPA Stimulator

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Vial of 1 mg/Flacon de 1 mg (Ref PP021B)  
Set of 6 vials of 1 mg/ Coffret de 6 flacons de 1 mg (Ref PP021K)

Ref. PP021B / PP021K

STORE AT  
2-8 °C

Lot :

Expiration :

CONSERVER A  
2-8 °C

### ANALYSIS CERTIFICATE

ANALYTICAL DATA	SPECIFICATIONS
1. <u>Protein Content (A280 nm)</u>  ..... mg/vial	> 0.8 mg/vial
3. <u>Reactivity in a tPA chromogenic assay</u>  Passed <input type="checkbox"/> Refused <input type="checkbox"/>	Passed

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DATE :

CONCLUSIONS :  Passed/Conforme     Refused/Non-conforme

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