HYPHEN BioMed

... a link between research and specialized applications for the biomedical diagnosis...

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STORE AT
2-8°C

For in vitro diagnostic use only

Diagnostic in vitro exclusivement

CONSERVER A
2-8°C

Composition:

Solution containing Calcium Chloride (0.25M) and Sodium azide (0.9 g/l) as preservative.

***CAUTION:** Sodium azide (NaN₃), may react with lead and copper plumbing to form highly explosive metal azides. In order to avoid this risk, flush with large volumes of water when discarding into a sink.

10 fold concentrated (M/4) Sodium Chloride solution

Presentation:

25 ml vial, 10 fold concentrated.

AR002A: 1 x 25 ml

AR002K: 3 x 25 ml

Use:

Calcium Chloride solution for coagulation assays, such as the Activated Partial Thromboplastin Time (APTT) or for the assays of factors of the intrinsic coagulation pathway.

10 fold concentrated solution. Dilute **1:10** with distilled water for obtaining a 0.025M (M/40) calcium chloride solution.

Note 1: Take care of stabilizing the solution at room temperature (18-25°C) for 30 minutes before use.

Note 2: Refer to the specific associated kit inserts for instructions concerning the tests to be performed.

For in vitro use only.

Stability:

In its original vial, unopened and kept at 2-8°C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8°C.

The 1:10 diluted solution, must be used within 48 hours.

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