
# HLA-B27 RealFast™ Assay

ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45

A-1120 Vienna, Austria

Phone: (+43-1) 8120156-0

[info@viennalab.com](mailto:info@viennalab.com)

[www.viennalab.com](http://www.viennalab.com)

## 1. Intended Use

The HLA-B27 RealFast™ Assay is a fast and accurate real-time PCR test for the detection of the HLA-B27 allele, a specific variant of the *human leukocyte antigen B (HLA-B)* gene that is strongly associated with seronegative spondyloarthropathies. The kit is used to confirm a suspected diagnosis of ankylosing spondylitis, reactive arthritis, juvenile rheumatoid arthritis or anterior uveitis. The qualitative assay discriminates the presence or absence of HLA-B27 in a human genomic DNA extract and detects all disease-relevant HLA-B27 subtypes. Reference sequence: HGVS: NG\_023187.1

## 2. Introduction

HLA-B molecules are cell surface proteins that play an important role in immunity. The frequency of the HLA-B27 variant is about 8% in the general Caucasian population, whereas it is found in over 90% of patients suffering from ankylosing spondylitis (AS). The disease, predominantly striking men between the age of 20 and 40 years, is characterized by inflammation of the sacroiliac joint and progressing stiffness of the spine. Due to the high correlation with the disease, HLA-B27 genotyping is suitable for the differential diagnosis of AS. Other subsets of spondyloarthritis are also associated with HLA-B27, although to a lesser degree.

## 3. Kit Contents

|                             |        |              |         |
|-----------------------------|--------|--------------|---------|
| RealFast™ 2x Genotyping Mix | 1 vial | □ white cap  | 1000 µl |
| HLA-B27 Assay Mix           | 1 vial | ■ purple cap | 550 µl  |
| HLA-B27 Positive Control    | 1 vial | ■ green cap  | 75 µl   |

The RealFast™ 2x Genotyping Mix comprises HotStart Taq DNA polymerase and dNTPs in an optimized buffer system.

The HLA-B27 Assay Mix consists of gene-specific primers and dual-labeled hydrolysis probes for HLA-B27 and a control gene.

A positive control for HLA-B27 is supplied with the kit.

The kit contains reagents for 100 reactions in a final volume of 20 µl each.

## 4. Storage and Stability

HLA-B27 RealFast™ Assay is shipped on cooling blocks. On arrival, store the kit at -20°C. Alternatively, store at 2 to 8°C for short-term use within one month. The kit withstands up to 20 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to intense light. If stored correctly, the kit will retain full activity until the expiration date indicated on the label.

## 5. Product Description

### 5.1. Principle of the Test

The test is based on the fluorogenic 5' nuclease assay, also known as TaqMan® assay. Each reaction contains gene-specific primer pairs which amplify a 184 bp fragment of the HLA-B27 gene and a control gene, the latter serving as PCR control. Further components are two dual-labeled, gene-specific hydrolysis probes which hybridize to the target sequence of the corresponding fragment. The proximity of the 5'-fluorescent reporter and 3'-quencher dye on intact probes prevents the reporter from fluorescing. During the extension phase of PCR the 5' – 3' exonuclease activity of the Taq DNA polymerase cleaves the 5'-fluorescent reporter from the hybridized probe. The physical separation of the fluorophore from the quencher dye generates a fluorescent signal in real-time, which is proportional to the accumulated PCR product. In samples positive for HLA-B27 both, the **FAM-labeled HLA-B27** probe as well as the **HEX-labeled PCR control** probe bind to the appropriate gene fragment. A strong fluorescence signal is detected in the FAM channel (520nm) and in the HEX channel (556nm). In samples negative for HLA-B27 only the HEX-labeled PCR control probe hybridizes to the complementary strand of the control gene fragment. A strong fluorescence signal is detected in the HEX channel and no or only a baseline signal in the FAM channel.

### 5.2. Real-time PCR Instrument Compatibility

The HLA-B27 RealFast™ Assay is validated for use with the AB 7500 Fast instrument.

The kit is compatible with various common real-time PCR instruments capable of recording FAM and HEX fluorescence:

- ✓ AB 7500 Fast (Applied Biosystems®)
- ✓ AB StepOne™ (Applied Biosystems®)
- ✓ CFX96™ (Bio-Rad)
- ✓ LightCycler® 480 (Roche)
- ✓ Mx3005P (Agilent Technologies)
- ✓ Rotor-Gene® 6000 (Qiagen)

» **Note:** RealFast™ Variant Detection QuickGuides for setting up and analyzing experiments on different types of instruments can be downloaded from [www.viennalab.com](http://www.viennalab.com).

When using AB StepOne™, set passive reference dye to "ROX ! «

The kit is supplied **without ROX**. For use with real-time PCR instruments requiring high ROX for normalization of data (e.g. Applied Biosystems® instruments StepOne™, 7300, 7900/7900HT), add ROX at a final concentration of 1 µM to the 2x Genotyping Mix.

### 5.3. Assay Performance Specifications

Determination of **sensitivity** was performed on 66 alleles testing positive for the HLA-B27 allele with a CE-marked reference kit. The HLA-B27 RealFast™ Assay determined all 66 alleles as positive, which equaled a true positive rate of 100%.

Determination of **specificity** was performed on 130 alleles testing negative for the HLA-B27 allele with a CE-marked reference kit. The HLA-B27 RealFast™ Assay determined all 130 alleles as negative, which equaled a true negative rate of 100%.

Limit of detection: 0.2 ng genomic DNA (per reaction)

Recommended DNA concentration: 2 to 20 ng/µl genomic DNA

## 6. Materials Required but not Supplied

Real-time PCR instrument with FAM (520 nm) and HEX (556 nm) filters, instrument-compatible reaction vessels, disposable powder-free gloves, vortexer, mini-centrifuge for 2.0 ml tubes, tube racks, set of calibrated micropipettes (0.5 – 1000 µl), sterile tips with aerosol-barrier filter, molecular grade water, DNA extraction system, freezer, biohazard waste container.

## 7. Experimental Protocol

### 7.1. DNA Extraction

DNA extraction reagents are **not supplied** with the kit.

DNA isolated from various specimens (e.g. whole peripheral blood, dried blood spots, buccal swabs or saliva) can be used. Ensure extracted DNA is suitable for amplification in terms of concentration, purity and integrity.

For accurate genotype calling, the DNA amount per reaction should be within the range of 10 to 100 ng for all samples.

### 7.2. PCR Controls

**Always** include a **No Template Control (NTC)** in each experiment to confirm absence of potential contamination. It is advisable to run the NTC (use PCR-grade water instead of DNA) in duplicate.

**Always** include the HLA-B27 **Positive Control** as positive reference signal for your unknown samples.

» **Note:** *The Positive Control is a potential source of contamination. Make sure to handle it carefully.* «

### 7.3. Preparation of HLA-B27 RealFast™ Master Mix:

Gently vortex and briefly centrifuge all solutions after thawing. Set up PCR at room temperature. Prepare sufficient **Master Mix** for all your reactions (N samples + positive control + negative controls) plus at least one additional reaction to compensate for pipetting inaccuracies:

| Component                   | per reaction | e.g. 24+1 reactions |
|-----------------------------|--------------|---------------------|
| RealFast™ 2x Genotyping Mix | 10 µl        | 250 µl              |
| HLA-B27 Assay Mix           | 5 µl         | 125 µl              |
| <b>Master Mix</b>           | <b>15 µl</b> | <b>375 µl</b>       |

Dispense **15 µl Master Mix** into each well. Add **5 µl purified DNA** or **Control** template to reach a final reaction volume of 20 µl.

To minimize risk of contamination, always pipette templates in the following order: first NTC, then samples, last positive control. Immediately close reaction vessels.

» **Note:** *Avoid creating bubbles in the final reaction mix and avoid touching the optical surface of the cap or sealing film without gloves. Both may interfere with fluorescence measurements. Centrifuge briefly if needed.* «

### 7.4. PCR Program

Program the real-time PCR instrument according to the manufacturer's instructions for quantitation experiments with two targets / reporter dyes. Place the samples into the thermal cycler and run the following program:

**AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480, Mx3005P**  
and other Peltier heating block-based instruments:

| Cycles | Temp        | Time   | Steps  |
|--------|-------------|--------|--|
| 1      | 95°C        | 3 min  | Initial denaturation   |
| 40     | 95°C        | 15 sec | Denaturation   |
|        | <b>60°C</b> | 1 min  | Annealing/Extension –<br><b>Data acquisition</b> on FAM<br>and HEX channel |

**Rotor-Gene® 6000:**

| Cycles | Temp                       | Time   | Steps  |
|--------|----------------------------|--------|--|
| 1      | 95°C                       | 3 min  | Initial denaturation   |
| 40     | 95°C                       | 15 sec | Denaturation   |
|        | 36-well rotor: <b>56°C</b> | 1 min  | Annealing/Extension –<br><b>Data acquisition</b> on<br>Green and Yellow<br>channel |
|        | 72-well rotor: <b>60°C</b> |        |  |

## 8. Data Analysis / Interpretation of Results


The presence or absence of the HLA-B27 allele is defined by whether there is a signal in the **FAM channel** or not. Successful PCR can be verified by an amplification of the control gene detected in the **HEX channel** (PCR control). Thus, genomic DNA samples positive for HLA-B27 as well as the HLA-B27 Positive Control show an amplification in both, the HEX and FAM channel. HLA-B27 negative samples show an amplification in the HEX channel only. Fluorescent levels and corresponding amplification curves are automatically displayed in amplification plots in the real-time PCR software.

| Sample Type              | Amplification in <b>FAM</b> channel (520 nm) | Amplification in <b>HEX</b> channel (556 nm) |
|--------------------------|--|--|
| HLA-B27 positive         | <b>YES</b>                                   | <b>YES</b>                                   |
| HLA-B27 negative         | NO   | <b>YES</b>                                   |
| HLA-B27 Positive Control | <b>YES</b>                                   | <b>YES</b>                                   |
| NTC                      | NO   | NO   |

Some instrument software may need manual threshold settings for accurate analysis.

Recommendations for Threshold Settings ( $C_0$ ):

Set threshold value for the FAM and for the HEX channel just above the background fluorescent signal generated by the No Template Control (NTC).

Amplification plot: **HLA-B27 positive**

## 9. Warnings and Precautions

- For *in vitro* diagnostics use only.
- Always use disposable powder-free gloves and wear suitable lab coat when handling specimens and reagents.
- Perform reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- Use pipettes dedicated for PCR setup only, use aerosol-guarded pipette tips.
- Use instrument-compatible reaction vessels with optically clear caps or sealers.
- Do not mix reagents from different lots.
- Do not use expired kits or kit components.

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**AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480, Mx3005P e outros instrumentos** de bloco de calor Peltier: