

CRYOcheck™

CorPac™

CRYOcheck CorPac consists of 30 individual vials containing human plasma with distinct PT, APTT and fibrinogen profiles spanning the normal and pathological range. Plasmas are buffered and rapidly frozen and are available in sets of 30 x 1.5 mL aliquots. Intended *For Research Use Only*.

DESCRIPTION	CATALOG #	FORMAT
CorPac	CCCP-15	30 x 1.5 mL

Expiry Dating: 2 years*



Normal Donor Set

CRYOcheck Normal Donor Set consists of platelet-poor plasmas prepared from carefully screened individual donors and is available in a set of 25 x 1.0 mL aliquots (50 unique donors available on request). Plasma is collected into trisodium citrate from healthy male and female donors between 18 and 66 years of age, processed, aliquoted and rapidly frozen. This results in a top quality, versatile product that truly represents “normal”. Intended *For Research Use Only*.

DESCRIPTION	CATALOG #	FORMAT
Normal Donor Set	CCNS-10	25 x 1.0 mL

Expiry Dating: 3 years*



Pooled Normal Plasma

CRYOcheck Pooled Normal Plasma is intended for use as a frozen alternative to in-house collected normal plasma pools and provides users with enhanced convenience and increased confidence in testing results. It may also be used as a normal control for prothrombin time (PT) and activated partial thromboplastin time (APPT). CRYOcheck Pooled Normal Plasma is comprised of platelet-poor plasma from 20 or more carefully screened male and female donors between 18 and 66 years of age. The plasma is buffered to enhance pH stability and frozen, resulting in a top quality, versatile product that truly represents “normal”.

DESCRIPTION	CATALOG #	FORMAT
Pooled Normal Plasma	CCN-10	80 x 1.0 mL
	CCN-15	80 x 1.5 mL
	CCN-40	81 x 4.0 mL
	CCN10-10	10 x 1.0 mL

Expiry Dating: 2 years*



- Platelet-poor (< 10,000 platelets/μL)
- Each lot number verified to be normal (factor assay value sheet available on request)
- 24-hour stability once thawed when maintained at 2 to 8 °C

CRYOcheck Features

- Donors are screened negative for lupus anticoagulants
- Source plasmas are screened negative for all FDA-required tests
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes

PrecisionBioLogic

www.precisionbiologic.com
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.

*From date of manufacture when stored at -40 to -80 °C