

CRYOcheck™
Factor V Deficient Plasma

Product Identification

Catalog Number: **FDP05-10**

Lot Number:

Expiry:

Assay	Value
PT	52.3 sec
APTT	157.7 sec
Fibrinogen	3.66 g/L
Factor II	100 %
Factor V	< 1.0 %
Factor V Antigen	< 0.01 U/mL
Factor VII	93 %
Factor VIII	102 %
Factor IX	117 %
Factor X	101 %
Factor XI	85 %
Factor XII	104 %

Results may vary depending on materials and methods used.

SOLUTION assured 

QD9 Rev.00 March 2007

CE

1.800.267.2796

info@precisionbiologic.com