

CRYOcheck™  
Factor II Deficient Plasma

**Product Identification**

Catalog Number: **FDP02-10**

Lot Number:

Expiry:

Assay	Value
PT	>240.0 sec
APTT	>240.0 sec
Fibrinogen	3.25 g/L
Factor II	< 1.0 %
Factor II Antigen	< 0.01 U/mL
Factor V	102 %
Factor VII	104 %
Factor VIII	109 %
Factor IX	138 %
Factor X	112 %
Factor XI	104 %
Factor XII	109 %

Results may vary depending on materials and methods used.

CE

SOLUTION assured

QD7 Rev.01 Oct 2016

1.800.267.2796

info@precisionbiologic.com