

CRYOcheck™
Factor IX Deficient Plasma

Product Identification

Catalog Number: **FDP09-10**
Lot Number:
Expiry:

| Assay | Value |
|-------------------|-------------|
| PT | 11.4 sec |
| APTT | 88.1 sec |
| Fibrinogen | 3.19 g/L |
| Factor II | 106 % |
| Factor V | 101 % |
| Factor VII | 112 % |
| Factor VIII | 94 % |
| Factor IX | < 1.0 % |
| Factor IX Antigen | < 0.01 U/mL |
| Factor X | 121 % |
| Factor XI | 104 % |
| Factor XII | 109 % |

Results may vary depending on materials and methods used.

CE

1.800.267.2796

info@precisionbiologic.com

SOLUTION *assured* ✓