

CRYOcheck™
Factor VII Deficient Plasma

Product Identification

Catalog Number: **FDP07-10**

Lot Number:

Expiry:

Assay	Value
PT	117.2 sec
APTT	32.1 sec
Fibrinogen	3.51 g/L
Factor II	97 %
Factor V	96 %
Factor VII	< 1.0 %
Factor VII Antigen	< 0.01 U/mL
Factor VIII	102 %
Factor IX	101 %
Factor X	102 %
Factor XI	95 %
Factor XII	102 %

Results may vary depending on materials and methods used.

CE

1.800.267.2796

info@precisionbiologic.com

SOLUTION assured 