

**CRYOcheck™**  
**Factor VIII Deficient Plasma**

**Product Identification**

Catalog Number: **FDP08-10**

Lot Number:

Expiry:

Assay	Value
PT	11.7 sec
APTT	81.5 sec
Fibrinogen	2.58 g/L
Factor II	108 %
Factor V	102 %
Factor VII	119 %
Factor VIII	< 1.0 %
Factor VIII Antigen	< 0.01 U/mL
Factor IX	125 %
Factor X	119 %
Factor XI	96 %
Factor XII	97 %

Results may vary depending on materials and methods used.

CE

1.800.267.2796

info@precisionbiologic.com

**SOLUTION** assured 