

CRYOcheck™
Factor X Deficient Plasma

Product Identification

Catalog Number: **FDP10-10**

Lot Number:

Expiry:

Assay	Value
PT	113.3 sec
APTT	99.0 sec
Fibrinogen	3.67 g/L
Factor II	107 %
Factor V	98 %
Factor VII	102 %
Factor VIII	101 %
Factor IX	105 %
Factor X	< 1.0 %
Factor X Antigen	< 0.01 U/mL
Factor XI	92 %
Factor XII	107 %

Results may vary depending on materials and methods used.

CE

1.800.267.2796

info@precisionbiologic.com

SOLUTION *assured* 