

CRYOcheck™
Factor XI Deficient Plasma

Product Identification

Catalog Number: **FDP11-10**

Lot Number:

Expiry:

Assay	Value
PT	11.9 sec
APTT	78.6 sec
Fibrinogen	3.40 g/L
Factor II	101 %
Factor V	104 %
Factor VII	109 %
Factor VIII	97 %
Factor IX	102 %
Factor X	114 %
Factor XI	< 1.0 %
Factor XI Antigen	< 0.01 U/mL
Factor XII	102 %

Results may vary depending on materials and methods used.

CE

1.800.267.2796

info@precisionbiologic.com

SOLUTION *assured* 