

CRYOcheck™ Lupus Positive Control

Product Identification

Catalog Number: **CCLP-10**
Lot Number:
Expiry:

Lupus Anticoagulant Assays	Reagent(s)	Test Results	Normal Range
APTT (Lupus Sensitive)	HemosIL APTT-SP®	99.0	26.1 – 43.1 sec
APTT (1:1 Normal Mix)	HemosIL APTT-SP®	72.1	26.1 – 43.1 sec
Kaolin Clotting Time (KCT)	Sigma	291	50 – 130 KCT units
Dilute Russell's Viper Venom Time (dRVVT)	CRYOcheck LA Check™ & LA Sure™	2.01	≤ 1.32 ratio
Hexagonal Lupus Anticoagulant (HLA)	Staclot® LA	46	0 – 8 sec
Tissue Thromboplastin Inhibition Time (TTI Ratio)	Thromborel S	1.5	0.5 – 1.2 ratio
Platelet Neutralization Procedure (PNP)	CRYOcheck Platelet Lysate	29.6	≤ 2.5 sec

The following results have been obtained using commercial ELISA assays, and are provided for additional information purposes only:

Analyte	Test Results
IgG Anticardiolipin Ab	Positive
IgM Anticardiolipin Ab	Negative
IgA Anticardiolipin Ab	Negative
IgG Anti-beta-2-glycoprotein 1	Positive
IgM Anti-beta-2-glycoprotein 1	Negative
IgA Anti-beta-2-glycoprotein 1	Negative
IgG Anti-phosphatidylserine	Positive
IgM Anti-phosphatidylserine	Negative

Results may vary depending on materials and methods used.

CE

SOLUTION *assured* ✓

1.800.267.2796

QD17 Rev.06 June 2012

info@precisionbiologic.com