

CRYOcheck™ Weak Lupus Positive Control

Product Identification

Catalog Number: **CCWLP-05**

Lot Number:

Expiry:

Lupus Anticoagulant Assays	Reagent(s)	Test Results	Normal Range
APTT (Lupus Sensitive)	HemosIL APTT-SP®	71.0	26.1 - 43.1 sec
APTT (1:1 Normal Mix)	HemosIL APTT-SP®	58.7	26.1 - 43.1 sec
Kaolin Clotting Time (KCT)	Sigma	200	50 - 130 KCT units
Dilute Russell's Viper Venom Time (dRVVT)	CRYOcheck LA Check™ & LA Sure™	1.59	≤ 1.32 ratio
Hexagonal Lupus Anticoagulant (HLA)	Staclot® LA	43	0 - 11 sec
Platelet Neutralization Procedure (PNP)	CRYOcheck Platelet Lysate	18.7	≤ 2.5 sec

The following results have been obtained using commercial ELISA assays, and are provided for additional information purposes only:

Analyte	Test Results
IgG Anticardiolipin Ab	Indeterminate
IgM Anticardiolipin Ab	Negative
IgA Anticardiolipin Ab	Negative
IgG Anti-beta-2-glycoprotein 1	Negative
IgM Anti-beta-2-glycoprotein 1	Negative
IgA Anti-beta-2-glycoprotein 1	Negative
IgG Anti-phosphatidylserine	Positive
IgM Anti-phosphatidylserine	Negative

Results may vary depending on materials and methods used.